

## G236p53DD/E1A Culture

Kazuhide Watanabe

6-21-2011

### Summary

G236p53DD/E1A was isolated from a mouse mammary tumor with MMTV-Wnt1/Py2<sup>flox/-</sup> genotype (mouse No. G236, B6/FVB/N mixed strain, established by KW). The isolated cells were cultured in 3D to form tumorsphere and then plated onto a plastic 2D dish, followed by the retroviral transduction of Dominant negative p53 (p53) and adenoviral E1A (pBABE-hygro p53 DD and pWZL hygro 12S E1A). Therefore, this line is hygromycin resistant. A bulk population of survived cells was established as a cell line. Pygo2 gene can be excised by exogenous Cre and Pygo2 protein was completely undetectable three days after Cre-adenovirus infection. A stable Wnt reporter line is also available. (lentiviral construct of TOP-flash 7TFP, puromycin resistant) (Fuerer and Nusse, 2010).

### Culture Medium preparation

1. In 1 liter beaker with a stir bar, add

800 ml milliQ water	
DMEM/F12 (1:1) powder (12500-062, Invitrogen)	1 bag
HEPES powder (BP310-500, Fisher)	5.95 g (25 mM final)
NaHCO <sub>3</sub> powder (S233-500, Fisher)	2.438 g
Pen-Strep aliquot (100x) (15140-122, Invitrogen)	10 mL

2. Adjust pH to 7.1-7.2 with 0.5M NaOH while stirring

3. Bring volume to 900 ml

4. Add the following:

Insulin (5 mg/ml in PBS) (I6634-250MG, Sigma)	1 ml	(5 µg/mL final)
EGF (100 µg/ml) (01-107, Millipore)	50 µL	(5 ng/mL final)
Hydrocortisone (4mg/ml in ethanol) (386698, Calbiochem)	250 µL	(1 µg/ml final)
Heat-inactivated FBS (varies)	100 ml	(10% final)

5. Mix well again by stirring

6. Filtered through 500 ml cup filter (0.22 µm, SCGP-T05-RE, Millipore) to TC bottles (200 ml/bottle)

7. Working media at 4 C; freeze unused media – can be stored in -20 °C for months before use

### **Routine culture (for 100-mm dish)**

Feed with fresh medium every 2-3 days. Usually split the cells every 3-5 days at 1:10 dilution.

**Note: these cells start dying rapidly after reaching confluency.**

1. wash twice with Ca<sup>2+</sup>-free PBS (5-10 mL).

2, add 2 mL of 0.25% Trypsin/EDTA and incubate for 3-5 min at 37°C.

3, Neutralize Trypsin by adding 5-mL of culture medium, then mix well, and centrifuge at 1000 rpm for 5 min.

4, aspirate supernatant and split the cells with 8-10 mL of fresh medium/dish.

### **Frozen stocks**

Resuspend 1-2 million cells in 1ml culture media with 10% DMSO, store in sealed Styrofoam at -80°C for at least overnight, then transfer to plastic box in liquid nitrogen tank.

### **Adenoviral infection (12-well plate)**

Day-1. Plate  $1 \times 10^5$  cells into a 12-well plate.

Day0. 1. Prepare  $5 \times 10^6$  TU (transduction unit, based on 293T cells) adenoviral solution by adding growth medium up to 200  $\mu$ L in total.

2. Aspirate medium from the cells and add adenoviral solution. Rock the plate to ensure that the entire surface area is covered.

3. Incubate the cells for 60 min in CO<sub>2</sub> incubator, rocking the plate every 15 min.

4. after 60 min, add 1 mL of growth medium.

Day1. replace medium or split cells.

Day2~ Confirm Pygo2 expression and perform experiments.

### **References**

Fuerer, C., and Nusse, R. (2010). Lentiviral vectors to probe and manipulate the Wnt signaling pathway. PloS one 5, e9370.